

Certificate

Certificate No.: BKR 23/10-12/15 Renewal decision dated: 03-10-2019

Expiry date: 03-12-2023

analysis methods
NF102 – Application to the food industry

The Company:

SOLABIA

41 rue Delizy 93698 Pantin Cedex - France

Is authorized to affix the NF VALIDATION mark on the alternative analysis method cited below, in accordance with the NF VALIDATION general rules and the certification rules NF102 - Validation of alternative analysis methods (Application to the food industry):

EASY STAPH®

Validated for the enumeration of coagulase-positive Staphylococci

Technical sheet reference's

EASY STAPH_V5 (fr)

This decision attests that the alternative analysis method has been assessed by AFNOR Certification and found to conform to the standards cited in page 2/2 and complementary requirements, as specified in the certification reference document. The **principal certified characteristics** are the "analytical performances" as defined in the associated validation study summarized report (relative trueness, accuracy profile...), available on the certification dedicated website http://nf-validation.org/en.

This certificate supersedes all previous certificates. This NF VALIDATION certificate, included 2 pages, is valid until **3**rd **December 2023**. It is subject to the results obtained upon regular controls carried out by AFNOR Certification. Appropriate decision is made by AFNOR Certification in accordance with the NF VALIDATION general rules and certification rules NF102 - Validation of alternative analysis methods (Application to the food industry).



Managing Director
Julien NIZRI

Issue dated 25/10/2021





Certificate No.: BKR 23/10-12/15 Renewal decision dated: 03-10-2019

Expiry date: 03-12-2023

The alternative analysis method:

EASY STAPH®

Validated for the enumeration of coagulase-positive Staphylococci

Manufactured by:

SOLABIA

Rue des 40 mines ZAC de Ther Allonne BP 10245 60002 Beauvais - France

Has been certified according to the reference documents and the application scope specified here after:

Validation protocol **NF EN ISO 16140-2** (September 2016): Microbiology of the food chain. Method validation - Part 2:

Protocol for the validation of alternative (proprietary) methods against a reference method.

Reference method(s) NF EN ISO 6888-2 (October 1999) and its amendment A1 (December 2003): Microbiology of food

and animal feeding stuffs. Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species) - Part 2: technique using rabbit plasma

fibrinogen agar.

Scope All human food products (by performing validation assays on a broad range of foods).

Restriction(s) None.

Warning None.

Other information Non applicable.

Please send any queries concerning the performances of the certified alternative method to AFNOR Certification.

You may download the validation study summarized report on http://nf-validation.afnor.org/en.

Issue dated 25/10/2021